

Our approach to safety

To create "safety" and "value"

As medical sophistication advances, the level of sophistication required in production, from research and development to manufacturing and shipping, is also increasing. YAMAKIN is strengthening relationships with external organizations and implementing world-class quality management, striving to provide safety and value.



Verification of Biological Safety



YAMAKIN established the Biological Science and Safety Laboratory in the Department of Oral and Maxillofacial Surgery, Kochi Medical School, Kochi University, Japan. We have been researching and analysing the safety risks of dental materials by the evaluation of various biological safety factors such as allergenicity, carcinogenicity, influence on our descendants, and so on.

YAMAKIN's Manufacturing Environment

All manufacturing processes for the TMR series are conducted in clean rooms and controlled extremely strictly.



Clean room in the YAMAKIN Kochi factory

Manufacturer

YAMAKIN CO., LTD.

1090-3 Otani, Kamibun, Kagami-cho, Konan-shi, Kochi, 781-5451 Japan https://www.yamakin-global.com Head Office: 3-7 Sanadayama-cho Tennoji-ku Osaka 543-0015, Japan Branch Office: Tokyo, Osaka, Sendai, Nagoya, Fukuoka, JAPAN Factory and R&D: Kochi, JAPAN P: +81-887-55-0281 F: +81-887-55-0053 E: contact@yamakin-gold.co.jp

C € 0123INTERNATIONAL20200928

INNOVATION Z Project No.1 Innovation by TMR series



Unique technology reduces the burden on the patients







Applicable from pulp capping to filling



Innovation by TMR Series






TMR series are developed and commercialized with ideals such as "no extracting", "minimal grinding" and "preserving teeth" based on the MI concept (minimal intervention).

TMR Series Technical Concept





Even though TMR-MTA cement Mielle, which will set with water, contains moisture, a compatible adherence surface can be assured, because TMR-AQUA BOND 0, which contains the phosphoric monomer M-TEG-P*, has amphiphilic property. Moreover, adherence is enabled even under moist conditions by the chemical interaction between zirconia, which is employed as a replacement for bismuth oxide, and the phosphate group of M-TEG-P.

Moreover, TMR-Z Fill 10. using YAMAKIN's filler technology, achieves both sustained fluoride release and high strength. Introduction the TMR series completes the total system from pulp capping to filling, with the long-term mission of "preserving teeth" in mind.



TMR-MTA cement Mielle

Improved X-ray radiography by 60%

Excellent usability with quick setting

- Easy mixing
- Japan made MTA cement features various advantages such as Bismuth Free.








TMR-AQUA BOND 0

Shorten the chair time.

Zero second bonding

- Stable adhesive strength in the wet condition
- Zero second demineralization reduces the burden on the patients.







TMR-Z Fill 10.

Achieves both sustained fluoride release property and strength.

- YAMAKIN's unique filler technology enables the both high strength and sustained fluoride release



TMR MTA cement Mielle


TMR-MTA cement Mielle Controlled medical device - Dental pulp capping materia

Excellent usability with quick setting "MTA cement" made in Japan through industry-academia collaboration




MTA cement with improved X-ray radiopacity

By increasing and optimizing the amount of zirconia, x-ray contrast or radiographic is improved by 60% compared to our previous product, making it easier to observe by x-ray and CT.



The percentage of calcium silicate (76vol%), which is the main ingredient, remains unchanged from the previous product. The percentage of alkaline and calcium ions releasing amount are also same. So, the basic performance of MTA cement, is maintained.





Bismuth Free

Zirconia is used as a radiographic agent for various biological cases, such as in dental materials and joint prosthesis.

Zirconia is a chemically stable substance so change of color is unlikely to happen.

		Exposure time of LED irradiator									
		0 second	10 seconds	30 seconds	90 seconds						
TMR MTA cement Mielle	White										
	Light Ivory										
Prototype co 20% of Bismu	-										

Note: Verified under oxygen insulated condition by glycerin.



Easy to Mix with a Small Amount of Water

At the moment the cement comes into contact with water, it is absorbed and becomes cement mud.

The bearing effect of the spherical silica particles improves the product's fluidity. This enables you to make an even paste easily with a small amount of water in a short time. Although we have increased the number of zirconia particles to improve the X-ray contrast, since the particles are the round shape, easy operation remains.






Image of bearing effect

High Compressive Strength

Since it sets quickly, the product achieves high compressive strength (about 90 MPa) from one day after application. The compressive strength reaches about 140 MPa after a week, matching the compressive strength of glass ionomer cement.



Moisture Ratio: 20 mass%

Quick Setting

Since it can be mixed with a small amount of water, the product sets quickly so that initial setting is completed in 15 to 30 minutes. Also, it sets sufficiently without any need to replenish the moisture with a moistened cotton pellet after application. If the paste gets dry and difficult to handle, the viscosity of the paste can be adjusted by adding moisture.

Moisture Ratio (Powder : Water)	State of Paste	Handling Time	Initial Setting Time		
20% (0.2g : 0.05g)	Standard	About 3 min.	15-30 min.		
25% (0.2g : 0.07g)	Soft	About 6 min.	30-40 min.		
30% (0.2g: 0.09g)	Very Soft	About 9 min.	40-60 min.		

Note 1: As indicated in the above chart, the initial setting time will be longer if the paste is softer. Remove excess moisture with a dry cotton pellet etc. after application.

Note 2: Note that in cases where the moisture ratio becomes 35% or more, remarkably the initial setting will become slow and the compressive strength will decrease.

Biocompatibility

New formation of dentine is confirmed, with odontoblast-like cells arrayed in a fence-like formation 1 week after direct pulp capping with TMR-MTA cement Mielle.





How to use TMR-MTA cement Mielle

Video clip for instructions on how to use MTA cement can be found on the web




TMR-MTA cement Mielle Product report



TMR-AQUA BOND 0

Controlled medical device - Dental Adhesive for Enamel and Dentine (Dental adhesive for ceramics) (Dental adhesive for metal)






▲∧M-TEG-P®



Advanced Adhesion and Speed

Moisture is necessary for demineralization and permeation of tooth substance in adhesion, but excess moisture is known to adversely affect adhesion. AQUA BOND 0 achieves consistent adhesion without sensitivity to moisture conditions, i.e., in both dry and moist conditions, with user-friendly application protocols due to the effects of one of the adhesive components, M-TEG-P.

An adhesive which prevents reduction of bonding strength caused by moisture (which cannot be removed clinically) has been developed.



▲∂M-TEG-P°


Fnamel Dentine



No Waiting Time for **Demineralization** (Tooth is demineralized quickly after application.)

One liquid, one step. A simplified adhesion step is employed. As time for demineralization is unnecessary, the burden on patients can be reduced.

AQUA BOND 0 can be used for various applications, such as children or elderly people, because of the shortened time required for operation in oral cavity.





No Separation After Taken Out

With precise control of ingredient blending and M-TEG-P, even after about 30 minutes under light-shielding conditions, there is no phase separation between lipophilic monomers and water; therefore, this enables the effective ingredients for adhesion to be applied uniformly.






M-TEG-P



Improves Bonding Strength of Composite

Filling minute roughnesses on the tooth surface caused by milling and creating smooth surface bonding layers enables the composite to adhere tightly. Furthermore, it enables a bonding layer (5 to 10 µm) that is as thin as possible and as uniform as possible, it fully demonstrates the physical strength of the composite.



Adheres to Various Materials

Thanks to M-TEG-P and other adhesive components, TMR-AQUA BOND 0 can bond a wide range of materials from tooth structure to zirconia, titanium and precious metals during CR treatment and intraoral repair. It also can adhere to porcelain, CAD/CAM produce crowns, and resin materials containing inorganic fillers during intraoral repair with using Multi Primer LIQUID as pretreatment.

	The substances in the table can be used										
	Tooth substance	Zirconia	Titanium	Gold alloy	Silver alloy	Gold-silver- palladium alloy	Resin (including inorganic material)	Porcelain material			
TMR-AQUA BOND 0	0	0	0	0	0	0	_	_			
Multi Primer LIQUID + TMR-AQUA BOND 0	_	_	_	_	_	_	0	0			

[Clinical case for using AQUA BOND 0 and Z Fill 10.]

Forming a Cavity

Form a cavity in accordance with normal

AQUA BOND 0 achieves consistent adhesion even in a moist condition.



2 Applying AQUA BOND 0

Apply AQUA BOND 0 to the whole surface of the cavity.

Air drying can be performed immediately after applying.

Light cure for 3 sec. with a LED light whose light intensity is 2,400 mw/cm².



3 Filling with Z Fill 10. Flow

Fill the cavity with the Flow type which has appropriate flowability. Many shades in the line-up enables the dental technician's image to be reproduced aesthetically.

Light cure for 4 sec. with a LED light whose light intensity is 2,400 mw/cm².

Restoration of cavity (Photos provided by: YAMAKITA DENTAL OFFICE (Kochi, Japan)



4 Filling with Z Fill 10. Universal

As the Universal type has good formability, delicate repair work is possible.


The A5 shade will be useful when there is less

Light cure for 4 sec. with a LED light whose light intensity is 2,400 mw/cm².



5 After the Treatment

Z Fill10. has both sustained fluoride release and strength, along with the aesthetic advantage of being like a natural tooth. This will improve the patient's QOL.







AQUA BOND CHALLENGE

Check it out the video clip showing the high adhesive strength of AQUA BOND 0.





TMR-AOUA BOND 0 Product report

An adhesive that has realized rapid treatment under moist conditions



Achieves Both Sustained Fluoride Release Property and Strength

TMR ZFIII 10.

Controlled medical device -Resin-based Dental Restorative Material







Achieves Both Sustained Fluoride Release Property and Strength

C.C.F.

Using this filler technology, not only is long-term sustained fluoride release achieved but also high strength and durability.

It has fluorine recharge characteristics, taking in fluorine contained in toothpaste and mouthwash and releasing it.





Dentistry-Polymer-based restorative materials.

Outstanding Usability

Universal is easily releasable from a spatula, but its viscosity allows it to stretch smoothly. In addition, as Flow and Low Flow have appropriate flowability, it is easy to fill or place them as desired.



Formability of Universal



Flowability of Flow (Left) Low Flow (Right)





TMR-ZFill10. CHAMELEON EFFECT






TMR-Z Fill 10. **Product report** Composite Resin: In Pursuit of Sustained Fluoride Release,

High Strength and Operability



Comprehensive Shade Line-Up



In combination with the filler technology, excellent light scattering property was achieved while maintaining the same transmittance as our conventional products.

A"chameleon effect" can be expected from this characteristic since the light scatters so that the surrounding colors become more natural. Furthermore, A5 and OA5 shades have been added to the line-up for Z Fill 10., and good aesthetics can be obtained even in cervical areas in elderly people or where the milling amount was small.

Shade line-up

Product name	Category	A1	A2	A3	A3.5	A4	A5	B1	B2	В3	C2	C3	D2	Others	QTY of shade	Contents
Z Fill 10. Universal	Dentine	•	•	•	•	•	•	•	•	•	•	•	•	BW, E, OW	17	3.8g (2mL)
	Opaque		•	•										DVV, E, OVV		
Z Fill 10. Flow	Dentine	•	•	•	•	•	•							DIAL F OIAL	14	2.6g (1.5mL)
	Opaque		•	•	•	•	•							BW, E, OW		
Z Fill 10. Low Flow	Dentine	•	•	•	•	•	•							DIA/ F OIA/	14	2.6g (1.5mL)
	Opaque		•	•	•	•	•							BW, E, OW		


BW: Bleaching White E: Enamel OW: Opaque White

YAMAKIN's Technology

Ceramics Cluster Filler

Ceramics cluster filler (C.C.F.) is YAMAKIN's original technology. C.C.F. consists of three kinds of filler whose uneven surfaces are intertwined together. By this "anchoring effect", it exhibits high strength. Moreover, the combination of the sustained fluoride release property and high strength, which is difficult to achieved in general by adding sustained fluoride release filler.

Note: more minute glass fillers than for C.C.F were used in Flow and Low Flow to give excellent flow control.



Ceramics Cluster Filler (1~20µm)

Spherical Nano-Filler (SiO₂: 20nm) Sustained Fluoride Release Filler

(Glass: Approx. 700 nm) Organic and inorganic composite filler

◀ Image of Z Fill 10. (Universal)



The concept of having both a sustained fluoride release property and high strength was first commercialized in HR blocks for CAD/CAM crowns: KZR-CAD HR BLOCK 2. Members of the development project were awarded the "Excellence Prize" for the product in the technological development category of the "Seventh Monodzukuri Nippon Grand Awards" by The Ministry of Economy, Trade and Industry (METI).

TMR ZFIII 10.

Direct Composite Resin

Controlled medical device -Resin-based Dental Restorative Material







Single Package

TMR-Z Fill 10. Universal 3.8g (2mL)

TMR-Z Fill 10. Flow 2.6g (1.5mL) Accessory, Needle Tip: 10 pcs.

TMR-Z Fill 10. Low Flow 2.6g (1.5mL) Accessory. Needle Tip: 10 pcs.

Set Package

Standard Pack (each 1 of A1, A2, A3, E)

· Universal • Flow

· Low Flow

A1x1, A2x1, A3x1, Enamelx1

2 pcs Pack (each 1 of A2, A3)

Universal

• Flow

· Low Flow

3 pcs of Same Shade Pack (Each 3pcs.of A2, A3, A3.5)

A2x1, A3x1

· Universal A2x3 • Flow

A3x3 · Low Flow A3.5x3

Accessory Needle Tip: 20 pcs.

Shade Card is included in each package.

Dental Adhesive

Controlled medical device

- Dental Adhesive for Enamel and Dentine (Dental adhesive for ceramics) (Dental adhesive for metal)





Single Package TMR-AQUA BOND 0 (5 mL)

TMR-AQUA BOND 0 Set

Set Package

• TMR-AQUA BOND 0 (5 mL) : 1 pc. • Disposable Applicator Brush: 50 pcs.

Disposable Plate: 25 pcs.

Set Package

2 pcs Pack of TMR-AQUA BOND 0

TMR-AQUA BOND 0 (5 mL) : 2 pcs.

Disposable Applicator Brush: 50 pcs. Disposable Plate: 50 pcs.



Packages and containers are subject to change without prior notice.

TMR MTA cement Mielle

Dental Pulp Capping Material

Controlled medical device - Dental pulp capping material





Single Package

TMR-MTA cement Mielle (White, Light Ivory)

• Microtube Type (0.2g): 3 pcs.

TMR-MTA cement Mielle (White, Light Ivory)

- Glass Container / Bottle Type (3g)
- · Accessory: Spoon (1 pc.)

TMR-MTA cement Mielle (White, Light Ivory)

- Glass Container / Bottle Type (10g)
- · Accessory: Spoon (1 pc.)



Multi Primer LIQUID

Bonding Material

 $\label{lem:controlled} \textbf{Controlled medical device} \ \ \textbf{-} \ \textbf{Bonding material for dental metal (dental ceramics, dental resin)}$



For intraoral repair of restorations

(metal, ceramics and resin)

For Repairing hybrid resin or porcelain prosthesis apply Multi Primer first, then use TMR-AQUA BOND 0 and TMR-Z Fill 10..

Single Package Multi Primer LIQUID (7mL)





Video clips





YAMAKIN CO., LTD.

