

Identification Number: 227AFBZX00042000
Controlled Medical Device – Ceramics materials

ZEO CE LIGHT ZR

[Contraindications and Prohibitions]

- Do not use this product on patients who have a history of allergy such as rash or irritation to crown restorations with this product or similar products.
- Do not use in patients with teeth grinding or oral parafunction.

[Product description]

[Intended purpose]

This product is used to produce restorations made from dental ceramics.

[Indications]

Fabrication of zirconia-ceramic restorations such as veneers, crowns, bridges, etc.

[Intended users]

Dental professionals only: dentists and dental technicians.

[Clinical Benefit]

Repair a patient's tooth by replacing the damaged part of the tooth as a zirconia-ceramic restoration.

[Shapes, Composition and Principles]

This product consists of the following items and contains the ingredients listed below:

Item	State	Ingredients
Liner	Paste	Glass powder in the system SiO ₂ , pigment, polyethylene glycol, etc.
Dentine	Powder	Glass powder in the system SiO ₂ , pigment, etc.
Enamel	Powder	Glass powder in the system SiO ₂ , pigment, etc.
Accent	Powder	Glass powder in the system SiO ₂ , pigment, etc.
Supplement	Powder	Glass powder in the system SiO ₂ , pigment, etc.
Stain	Paste	Glass powder in the system SiO ₂ , pigment, propylene glycol, etc.

[Intended use and effect / efficacy]

This product is used to produce restorations made from dental ceramics.

[Product Specification]

Physical Property (Test method: ISO6872:2015)

Item	Requirement
Flexural strength	More than 50 MPa
Chemical solubility	Max. 100 µg/cm ²
Liner thermal expansion	Shall not deviated by more than 0.5 × 10 ⁻⁶ /K from the stated value Stated value: Firing 2 times, 4 times Liner: 10.0 Dentine, Enamel: 9.6
Glass transition temperature	Shall not deviated by more than 20°C. Degree from the stated value. Start Value: Firing 2 times, 4 times Liner: 610°C Dentine, Enamel: 590°C

[Firing Schedule]

	Dry Time (min.)	Start Temp. (°C)	Heat Rate (°C/min.)	Firing Temp. (°C)	Hold Time	Reducing Pressure (kPa)
Liner	5-10	450-500	60	920-940	1	97
Dentine	6-9	550-600	60	900-920	1	97
Self-Glaze	4	550-600	60	880-900	1	Air
Supplement	6	550-600	60	800-850	1	97
Stain	4	550-600	45-60	850-870	1	97

(NOTE)

- In case where Stain is used on the surface of Zirconia, Heat Rate should be set 45°C/min. In case where Stain is used for Ceramics, it should be set 60°C/min.
- Conditions, such as the above schedule is a guideline, not precise data. Conditions of firing are different depending on kinds or shapes of furnaces. Confirm the temperature and time by test firing before clinical case.

[Operating Instructions]

1. Liner Application and Firing

- In case where using Liner, paste type, agitate Liner in the container, scoop requisite amount on the pallet and then mix it until it reaches a uniform consistency by a spatula.
- Apply Liner thinly and condense lightly then fire as above mentioned schedule.

2. Dentine, Enamel and Accent Application and Firing

Mix each Dentine, Enamel, Accent with distilled water, and build up to the crown form. Fire as above mentioned schedule after condensation and water absorption.

3. Correction to shape and Glaze Firing

- Correct the shape using a carborundum, diamond point or white point.
- After correction to shape, clean by an ultrasonic cleaner with diluted water or alcohol. Then fire for glaze.
- If necessary, build up by applying Stain or Supplement and fire.

[Technical Directions and Points for Attention]

- Do not use this product in combination with other products.
- Do not grind or polish this product in the mouth.
- This product should be used on zirconia restoration.

[Precautions]

1. Precautions

- Stop using this product and seek medical attention if any sign of allergy, such as rash, eczema, redness, ulcer, swelling, itching, numbness and etc. appears in patients.
- Operators who have a history of allergy such as irritation or rash to this product or organic solvent must wear plastic gloves to avoid direct skin contact. Operators must discontinue using and seek medical attention if any sign of allergy, such as rash, eczema, redness, ulcer, swelling, itching, numbness and etc.
- In case of eye contact, immediately flush with running water and consult a physician.

2. Important Basic Cautions

- When grinding and polishing this product, use exhaust fan and anti-dust mask approved by your local public safety agency in order to avoid inhaling dust. Wear safety goggles to protect the eyes.
- Provide adequate ventilation during operation.

[Handling and Storage]

Keep out of reach of personnel other than dental staff.

[Expiry Date]

- This product must be used before the expiry date printed on the package.
- The expiry date printed on the package is based on our authentication.
- The expiry date printed on the package refers to use expiry date.

* (ex.: YYYY-MM means the last day of YYYY (year)/MM (month) as the expiry date)

[ZEO CE LIGHT ZR: Package]

Item	Shade	State
Liner	L-A(1,2,3,3.5,4), L-B(1,2,3,4), L-C(1,2,3,4), L-D(2,3,4), L-G(1,2,3,4,5,Or), L-W(1,2,3,4)	Paste
Dentine	D-A(1,2,3,3.5,4), D-B(1,2,3,4), D-C(1,2,3,4), D-D(2,3,4), D-W(1,2,3,4)	Powder
Enamel	E-(1,2,3,4), OE(1,2,3,4,Amber)	Powder
Accent	Opal 1, Opal 2, Snow White, White, Vanilla, Milky, Yellow Ocher, Mamelon Vanilla, Mamelon Orange, Mamelon Yellow, Occlusal Brown, Occlusal Orange, Occlusal Enamel, T-Natural, T-Glass, T-Clear, T-Opacity, T-Buffer, T-Amber, T-Salmon, T-Gray, T-Violet, T-Blue, T-Baby Pink, T-Pink, G(1,2,3,4,5,6,7,Or)	Powder
Supplement	SUP-(D,E,T)	Powder
Stain	Glaze, White, Milky, Light Orange, Orange, Brown, Dark Brown, Gray, Violet, Blue, Pink, Yellow, Black, Indigo, Dark Pink, Amber, Khaki, Buff, Fluorescent, Opal Glaze, Shade A, Shade B, Shade C, Shade D	Paste

 YAMAKIN CO., LTD. 1090-3 Otani, Kamibun, Kagami-cho, Konan-shi, Kochi, 781-8451, Japan E-Mail : contact@yamakin-gold.co.jp	CE 0123  MT Promedt Consulting GmbH Ernst-Heckel-Straße 7, 66386 St. Ingbert, Germany Tel : +49 6894 581020 Fax : +49 6894 581021 E-Mail : info@mt-procons.com
---	--