

Save Your Time with



Luna-Wing
Hybrid Composite Resin for C&B

&

LED CURE Master
LED Light Curing Machine




50min.
Porcelain
Fused to Metal

Fabricating Time*

*Fabricating Time is estimated by calculating only firing time and curing time for porcelain and Luna-Wing respectively.

2min. 20sec.
Luna-Wing

Approximate Firing/Curing Time (tested at YAMAKIN)				
	Opaque	Dentine	Enamel	Total
PFM	30min.*	10min.*	10min.*	50min.*
 Luna-Wing	40sec.**	10sec.**	90sec.**	2min 20sec.**

*The time varies depending on furnaces **Cured with LED CURE Master

For Higher Properties

TWiNY

High-Class Indirect Composite Resin



For more information, go to the website at
▶ www.yamakin-global.com/twiny




Luna-Wing


Hybrid Composite Resin for C&B

- ✓ Easy to Handle!
- ✓ Extremely Low Shrinkage!
- ✓ Good Shape Preservation!
- ✓ Nano-Technology!
- ✓ High Biological Safety!

Basic Layering Steps for Luna-Wing



1 Apply Primer Paste & Invisible Opaque + Light Cure



2 Apply Opaque+ Light Cure



3 Apply Dentine + Light Cure




4 Apply Enamel + Light Cure



5 Polish & Glaze



6 Finish



Extensive Color Variation

Basic Shade Color (A0 - D4)

A0	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
MO																
OA0	OA1	OA2	OA3	OA3.5	OA4	OB1	OB2	OB3	OB4	OC1	OC2	OC3	OC4	OD2	OD3	OD4
CA1 CA2 CB1 CB2 CC1 CC2 CD1 CD2																
ODA0	ODA1	ODA2	ODA3	ODA3.5	ODA4											
DA0	DA1	DA2	DA3	DA3.5	DA4	DB1	DB2	DB3	DB4	DC1	DC2	DC3	DC4	DD2	DD3	DD4
E0	E2	E3	E4	E1	E2	E3	E2	E3	E4	E2	E3	E4	E2	E3	E4	E3

Effect Line-Up 9 colors



+ Stain Line-Up 21 colors

