

Controlled Medical Device – Dental glazing material (Coloring material for polymer-based crown, Bonding Material for Dental Resin)

## Nu:le Coat

Characterization material for resin

Pre-treatment for PEEK

# Nu:le Coat

More easily, more satisfactory.



Manufacturer

**YAMAKIN CO., LTD.**  
1090-3 Otani, Kamibun, Kagami-cho,  
Konan-shi, Kochi, 781-5451 Japan  
<https://www.yamakin-global.com>

Head Office: 1090-3 Kamibun, Kagami-cho, Konan-shi, Kochi, 781-5451, Japan  
Branch Office: Osaka, Tokyo, Nagoya, Fukuoka, Sendai, JAPAN  
Factory and R&D: Kochi, JAPAN  
P: +81-887-55-0281  
E: [contact@yamakin-gold.co.jp](mailto:contact@yamakin-gold.co.jp)

**CE 0123**

INTERNATIONAL20231213



For Laboratory  
&  
Clinic Use



Low Viscosity for  
Easy Application



Transparent Glossy Coat  
Can be overglazed



Excellent  
Durability



For Flexible Color  
Adjustment



Bonding for  
PEEK frame


# Nu:le Coat



## Easy Application and Solid Coating


### Pursuing Transparency

The higher surface hardness remained the same and reduced the amount of photopolymerization initiator that causes yellowness to the extreme. Nu:le Coat Liquid Clear (hereafter Clear) is clear even immediately after light curing, so it is difficult to affect the original color tone of dental restoration even when applied overglaze. Therefore, Clear can be used easily for whitening shades.




### Small Curing Strain

The combination of flexible polyfunctional monomers and tough cross-linkable monomers reduced curing strain. Therefore, even if the applied layer become thicker because of overglaze, it is less likely to occur cracks or floated surface.




#### Crack Evaluation Test

Dropped the liquid directly into the pit and fissure of crowns and evaluated by microscopy for the presence of cracks and floated surface after light curing.



Thanks to the small curing strain, the inhibition of cracking and floated surface can be expected even when a lot of liquid has been applied to the pit and fissure or the junction.

#### Evaluation of Curing Strain




Curing deformation after light curing (strain amount) is greatly reduced comparing the thickness before light curing (0.5 mm).

## Excellent Durability

In the evaluation of toothbrush abrasion test (50000 times) <sup>\*1)</sup>, it is equal to or more than that of polishing finish, and excellent abrasion resistance has been confirmed.

It is reported that plaque adherence is rapidly increase when the surface roughness exceeds 0.2  $\mu\text{m}$  <sup>\*2)</sup>. However, the surface roughness of Nu:le Coat applied surface is less than 0.1  $\mu\text{m}$  even after 50,000 times of toothbrush abrasion test.

<sup>\*1)</sup> ISO14569-1, the ratio of toothpaste and water is 1:2, Load weight 2.0N  
Brush stroke 50,000 times is equivalent to about 7 years if it is assumed to brush teeth twice a day and brush 10 times per tooth.  
<sup>\*2)</sup> Bollen. M. Et al. : Dental Mater, 13 (4) 258~269, 1997.



## Low Viscosity and Easy Application

Nu:le Coat Liquid was pursued an optimal viscosity that is easy to apply. Thanks to low viscosity design, hard to remain traces of a brush and easy to apply thinly and evenly. Thin hard-coated layers of about 5  $\mu\text{m}$  can be created, allowing coating and characterization without losing the surface structure of dental restorations.

And please use Gel type when thicker layers are necessary.

#### Thickness of 1 layer


- Liquid : 5-15  $\mu\text{m}$
- Gel : 20-50  $\mu\text{m}$



Easy glazing with easy steps. There is no loss or chipping concern during polishing.

## Shorten the working time with Nu:le Coat


Finish with polishing requires craftsmanship so polishing is a very valuable task. However, depending on the level of proficiency, the finish may vary. On the other hand, Nu:le Coat is easy to apply and provide glaze on the surface in a short period of time regardless of proficiency. The working time becomes 1/3 to achieve glaze on the surface.



## Get More Aesthetic Qualities With Color Type


### Wide Range of Color Line-up for Crown Characterizing

White, Gray, Brown, Yellow, Orange and Blue lineups allow for free characterization of crowns.



### Shade plus

There are A plus, B plus, C plus, and D plus lineups that allow fine adjustment of the shade.




A2 Crown Applied 1 layer of A plus.



B2 Crown Applied 1 layer of B plus.




C2 Crown Applied 1 layer of C plus.



D2 Crown Applied 1 layer of D plus.

### Gum Shades

There are also Gum shade line-up (Red, Pink and Black) for denture application or Gum expression.




Characterized resin frame by NuLe Coat Liquid Red, Pink, Black and Orange

## Plus One Application

### Pre-treatment of PEEK before applying resin materials.

Super engineering plastics PEAK-based (PEEK, PEKK, etc.) have attracted attention in recent years for their ability to mitigate impact on bone in implant cases. By the combination of Clear and Gel, composite resin (e.g., Luna-Wing) and hybrid composite resin (e.g., TWiNY) can be built up on PEEK frame.

#### Shear Bond Test (PEEK and Opaque Resin)



Applied Clear and Gel (at 1:1 ratio) after sandblasting

Completion after TWiNY application




## For Both Dental Laboratory and Chairside Use

Nu:le Coat can be used not only for the finishing of resin-based dental restorations but also for regular maintenance of resin-based dental restorations that have lost their glaze.

### Procedure of Applying Nu:le Coat

#### For Laboratory Use



After shape modification    Roughen the surface (Alumina Sandblast etc.)    Application of Nu:le Coat (characterizing and glazing)    Light Curing (60 sec. by LED CURE Master)    Completion


#### For Chairside Use



Filling composite    After light irradiation, shape modification and roughening the surface.    Application of Nu:le Coat (glazing)    Light irradiation    Completion

\*This product is for resin restoration only. This product cannot be used for tooth substances.

**Point** Gel is recommended to use for the parts where need thickness.



For Color types, shake the bottle well before use as the colouring material may settle. Please avoid using this product for occlusal parts where pressure is applied.

Nu:le Coat  
Movie



## Line-up

Nu:le Coat

Nu:le Coat

Controlled Medical Device – Dental glazing material  
(Coloring material for polymer-based crown, Bonding Material for Dental Resin)



Single Package

Nu:le Coat Liquid Clear (6ml)  
Nu:le Coat Liquid (4ml) 15 colors  
(White, Gray, Brown, Yellow, Orange, Blue, Red, Pink, Gum, Dark Brown HV, Black, A plus, B plus, C plus, D plus)



Single Package

Nu:le Coat Gel (2ml)

Set Package

#### Intro Set

- Nu:le Coat Liquid Clear (6ml) : 1 pc
- Nu:le Coat Gel (2ml) : 1 pc
- Round Brush : 1 pc
- Flat Brush : 1 pc
- Disposal Plate : 10 pcs

#### Regular Set

- Nu:le Coat Liquid Clear (6ml) : 1 pc
- Nu:le Coat Gel (2ml) : 1 pc
- Nu:le Coat Liquid (4ml) : White, Gray, Blue, Brown (1 pc each)
- Round Brush : 1 pc
- Flat Brush : 1 pc
- Disposal Plate : 20 pcs

#### Full Set

- Nu:le Coat Liquid Clear (6ml) : 1 pc
- Nu:le Coat Gel (2ml) : 1 pc
- Nu:le Coat Liquid (4ml) : White, Gray, Blue, Brown, Yellow, Orange, A Plus, B Plus, C Plus, D Plus (1 pc each)
- Round Brush : 1 pc
- Flat Brush : 1 pc
- Disposal Plate : 30 pcs

#### Set Contents

| Type        | Gel   |       | Liquid |      |       |        |        |      |     |      |     |               |       |        |        |        |        |  |
|-------------|-------|-------|--------|------|-------|--------|--------|------|-----|------|-----|---------------|-------|--------|--------|--------|--------|--|
|             | Clear | Clear | White  | Gray | Brown | Yellow | Orange | Blue | Red | Pink | Gum | Dark Brown HV | Black | A Plus | B Plus | C Plus | D Plus |  |
| Intro Set   | ●     | ●     |        |      |       |        |        |      |     |      |     |               |       |        |        |        |        |  |
| Regular Set | ●     | ●     | ●      | ●    | ●     |        |        | ●    |     |      |     |               |       |        |        |        |        |  |
| Full Set    | ●     | ●     | ●      | ●    | ●     | ●      | ●      | ●    |     |      |     |               |       | ●      | ●      | ●      | ●      |  |

#### Light Curing Time

| Usage           | Equipment*1         | Time   |              |
|-----------------|---------------------|--|--------------|
|                 |                     | Pre-curing   | Final curing |
| Laboratory Use  | LED CURE Master     | Pre-curing   | 30 sec.      |
|                 |                     | Final curing                                       | 60 sec.      |
|                 |                     | Pre-treatment for PEEK                             | 30 sec.      |
| Chairside Use*2 | Dental Curing Light | Pre-curing STD 10 sec.<br>Final curing STD 20 sec. |              |

\*1) When using a photopolymerization device other than those listed above, refer to IFU for each system and check that the effective wavelength (e.g. purple LED: peak output wavelength 400 ± 20 nm) is suitable for the material before use. For the effective wavelength, please contact the distributor of the photopolymerization device.

\*2) For intraoral use and large application areas, multiple light irradiations should be performed.

#### Related Products



Luna-Wing  
Hybrid Composite Resin for C&B



TWiNY  
Indirect composite resin



LED CURE Master  
LED light cure apparatus